GAS **S** **TEST**® **FOUR VECT**

Ref. GT010204
Ref. GT010205
Ref. GT010210
Ref. GT010220

Sólo para uso veterinario / For veterinary use only.

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E Teske, F van Knapen, EGM Beijer, RJ Slappendel: " Risk of Infection with Leishmania spp. in the Canine Population in the Netherlands" in Acta Veterinaria Scandinavica 2002, 43:195-201.

Lindsay DS, Zajac AM, Barr SC: Leishmaniasis in American Foxhounds: An Emerging Zoonosis? Compend Cont Educ Pract Vet 24:304-312, 2002.

Torsten J. Naucke: „Leishmaniose, eine Tropenkrankheit und deren Vektoren (Diptera, Psychodidae, Phlebotominae) in Mitteleuropa“ in Denisia 0006, 2002 (163-178).

Al-Adhami B., Scandrett W.B., Lobanov V.A. & Gajadhar A.A. (2011): Serological cross-reactivity between Anaplasma marginale and an Ehrlichia species in naturally and

experimentally infected cattle. *Journal of Veterinary Diagnostic Investigation* 23, 1181-1188.

Alleman, A. R. and Wamsley, H. L. :An update on anaplasmosis in dogs. *Vet Med* 2008, April: 212-220.

CDC: Ehrlichiosis and Anaplasmosis: 2008 Case Definition. US Department and Health and Human Services, CDC, National Notifiable Diseases Surveillance System (NNDSS): Atlanta, GA.

Dujardin J-C, Campino L, Carmen C, Dedet, J-P, Gradoni L, Soteriadou K, Mazeris A, Ozbek Y and Boela M. 2008. spread of vector-borne diseases and neglect of leishmaniasis, Europe. *Emerging Infectious Diseases*. *Emerging Infectious Diseases*. 14, no. 7.

Eckert, Johannes, Friedhoff, Karl Theodor, Zahner, Horst, Deplazes, Peter. Textbook of parasitology for veterinary medicine. 2nd, completely revised edition. Enke-Verlag, 2008 European Scientific Counsel Companion Animal Parasites: Control of vector-borne diseases in dogs and cats.

Maggi, R. G.; Birkenheuer, A. J.; Hegarty, B.C.; Bradley, J. M.; Levy, M. G.; Breitschwerdt, E. B. : Comparison of serological and molecular panels for diagnosis of vector-borne diseases in dogs. *Parasite Vectors*. 2014;7:127.

Pantchev, N.: Tick-borne travel infections in dogs: ehrlichiosis and babesiosis. In: *Veterinary Review*. Vol. 4, 2012, pp. 162-170.

Ranjbar-Bahadori, S. H.; Veshgini, A.; Shirani, D.; Eslami, A.; Mohieddin, H.; Shemshadi, B.; Masooleh, R. :Epidemiological aspects of canine dirofilariasis in the north of Iran. *Iran J Parasitol*. 2011 Mar;6(1):73-80.

Sainz, Á.; Roura, X.; Miró, G.; Estrada-Peña, A.; Kohn, B.; Harrus, S.; Solano-Gallego, L.: Guideline for veterinary practitioners on canine ehrlichiosis and anaplasmosis in Europe. *Parasite Vectors* 2015 Feb 4;8:75. doi: 10.1186/s13071-015-0649-0.

Self, S.; Liu, Y.; Nordone, S. K.; Yabsley, M. J.; Walden, H. S.; Lund, R. B.; Bowman, D. D.; Carpenter, C.; McMahan, C. S.; Gettings, J. R.: Canine vector-borne disease: mapping and the accuracy of forecasting using big data from the veterinary community. *Anim Health Res Rev*. 2019 Jun;20(1):47-60. doi: 10.1017/S1466252319000045.

Vieira, A. L.; João Vieira, M.; Oliveira, J. M.; Simões, A. R.; Diez-Baños, P.; Gestal, J. : Prevalence of canine heartworm (*Dirofilaria immitis*) disease in dogs of central Portugal. *Parasite*. 2014; 21:5. doi: 10.1051/parasite/2014003. epub 2014 Feb 19.

1. Intended Use

The GaSS-Test FOuR VECt test is used to detect vector diseases such as Leishmania, Ehrlichia, Anaplasma, detecting antibodies, and Dirofilaria, detecting antigens, in samples obtained from dog blood. This test is designed to aid in the diagnosis of the four diseases as with an average sensitivity of 97% and specificity of 99.99% it is highly likely to detect Leishmania infantum (Ac), Ehrlichia canis (Ac), Anaplasma phagocytophilum (Ac) and Dirofilaria immitis (Ag).

2. Introduction and clinical significance

Vector-borne diseases in dogs are becoming increasingly important due to various factors. Their occurrence is no longer limited to individual countries, but is an increasingly global problem (European Scientific Council for Companion Animal Parasites (ESCCAP).

Due to climate change arthropods (ticks, diptera, sandflies etc.) are adapting to changing temperatures, so that ticks, for example, are no longer only a problem in summer, but are also active in winter.

In addition, mobility is increasing worldwide: dogs are being taken on trips to other countries and also more and more dogs are finding a new home in another country or on another continent (Self et al., 2019).

Vector-borne diseases such as Leishmania, Ehrlichia, Anaplasma and Dirofilaria can severely damage the

body and even cause death. Therefore, reliable detection methods should be used in suspected cases and the results should always be checked with a veterinary reference laboratory.

The GaSS-Test FOuR VECt has been developed for this purpose, it is a rapid test based on the principle of the lateral flow system. The test detects the following relevant pathogens in the blood, serum or plasma of dogs: Leishmania infantum (Ab), Ehrlichia canis (Ab), Anaplasma phagocytophilum (Ab) and Dirofilaria immitis (Ag).

All pathogens are known as Mediterranean diseases, but by now the vectors - mosquitoes (including Culex, Aedes, Anopheles and ticks including Ixodes ricinus) - are also native to other areas (Maggi et al., 2014).

2.1 Pathogens

Leishmania: Leishmania infantum is a serious veterinary problem throughout the Mediterranean basin, where dogs are the most important reservoirs of the parasite. More than 500 species of sandflies are known, but only about 30 have been found to transmit leishmaniasis. Of these, only female sandflies transmit the parasite. These females need blood to develop their eggs and become infected with Leishmania parasites when they suck blood from an infected person or animal. Infected dogs can themselves develop leishmaniasis and are a potential source of contamination for other dogs

and humans. The parasites develop over a period of 4 to 25 days in the phlebotomine sandfly. When the infective female sandfly feeds on fresh blood, she inoculates the parasite into the person or animal and the transmission cycle is complete.

During canine leishmaniasis (CL) due to *Leishmania infantum*, high levels of antibody production are associated with the presence of various clinical signs, due to the deposition of soluble immune complexes in organs and tissues. Antibodies can usually be detected in these clinically affected animals. *Leishmania* antibodies can be successfully detected up to weeks or months after infection. Detection of specific antibodies in serum samples is widely used for the diagnosis of canine leishmaniasis. Asymptomatically infected animals often have marginal or very low antibody titers and are difficult to detect.

Ehrlichia canis: Dogs infected with *Ehrlichia* develop canine monocytic ehrlichiosis. The occurrence of infections is related to the spread of the vector *R. sanguineus*. For the diagnosis of *Ehrlichia* infections, a combination of anamnesis, to assess for possible tick infestation, evaluation of clinical signs, serology and/or PCR is available. Serologically, antibodies can be detected with the FOuR VECt test using *E. canis* specific antigens.

Anaplasma phagocytophilum: In Europe, *A. phagocytophilum* and *A. platys* infections occur in dogs. For detection, a mixture of *A. phagocytophilum* and *A. platys* antigens and synthetic

peptides of the immunodominant p44 surface protein are applied to the test strip. The use of purified synthetic peptides allows cross-reactivity of *A. platys* pathogens, and *A. phagocytophilum*. For diagnostic purposes, it is important to note that some studies have indicated that *A. phagocytophilum* strains are more virulent in the United States than strains circulating in Europe. Therefore, symptomatic infection in European patients appears to be rarer and less severe than in North America. Therefore, clinical investigation of patients with positive antibody detection is of great importance.

Dirofilaria immitis: *D. immitis* infection can cause severe and fatal disease in dogs. Adult heartworms can survive for 5-7 years in dogs. The antigen test is used to detect circulating *D. immitis* antigens in the blood; adult female heartworms are detected.

3. Principle of use

The GaSS-Test *Leishmania* is a lateral flow immunoassay. The test uses specific antigens that capture *Leishmania* antibodies from the sample allowing the agent to be visualized. The GaSS-Test *Leishmania* is a highly sensitive immunoassay and consists of a test strip housed in a convenient test cassette.

The GaSS-Test *Leishmania* test strips contain a sample pad, a conjugate pad, a reaction membrane and an absorption pad. After application of the sample to the test strip, the

Leishmania antibodies present in the sample material react with the labelled antigens, forming a complex on the conjugate pad. The liquid then migrates along the strip by capillary action. In the test line region, the complexes formed are captured by the immobilised antibodies, resulting in the formation of a red line in the test region indicating a positive result. If the sample does not contain Leishmania antibodies, the line in the test region will not appear. In addition, a red line should form in the control region, regardless of the presence of leishmania antibodies in the sample. This line serves as an internal control and indicates that an adequate volume of sample liquid has been added to the test strip and that the membrane is sufficiently soaked.

4. Materials provided

4.1 4 test kit, Art. Nr. GT010204

- 4 GaSS-Test FOuR VECt Tests
- 4 EDTA K3 tubes 0.5 ml
- 4 Disposable dropper for sample collection
- 2 Diluent (buffer) 4 ml
- 1 instructions for use

4.2 5 test kit, Art. Nr. GT010205

- 5 GaSS-Test FOuR VECt Tests
- 5 EDTA K3 tubes 0.5 ml
- 5 Disposable dropper for sample collection
- 2 Diluent (buffer) 4 ml
- 1 instructions for use

4.3 10 test kit, Art. Nr. GT010210

- 10 GaSS-Test FOuR VECt Tests
- 10 EDTA K3 tubes 0.5 ml
- 10 Disposable droppers for sample collection
- 2 Diluents (buffer) 4 ml
- 1 instructions for use

4.4 20 test kit, Art. Nr. GT010220

- 20 GaSS-Test FOuR VECt Tests
- 20 EDTA K3 tubes 0.5 ml
- 20 Disposable dropper for sample collection
- 4 Diluent (buffer) 4 ml
- 1 instructions for use

5. Materials required but not provided

Timer.

6. Storage and stability

The GaSS-Test FOuR VECt should be stored at room temperature (2-30 °C). The test remains stable until the expiry date printed on the package. Keep the cassette in its packaging until use. Do not freeze the tests or subject them to temperatures above 30 °C.

7. Precautions

- Use a new test cassette for each test.
- Single use only, do not reuse the test.
- For veterinary use only.

- Use only the materials provided with the GaSS-Test kit.
- Use the test within 60 minutes after opening the pouch.
- The test cassette must lie flat on a smooth surface.
- Note the amount of sample material used. Incorrect number of drops or drops that are too small may lead to false results.
- Consider the test results as invalid after the reading time.
- Do not use the kit when its pouch and the sealing are damaged or the expiration date has passed.
- Dispose of all contaminated materials properly according to standard regulations for handling potentially infectious materials and chemical reagents during testing. Protective clothing is recommended. Disinfect the work area after testing. Dispose of sample materials according to local regulations.

8. Test procedure

Serum and plasma samples


Recommended sample material is a freshly collected serum or plasma to achieve the highest detection sensitivity.

Separate serum or plasma from whole blood as quickly as possible, the sample should be clear and not haemolysed or lipaemic.

Test procedure

Open the pouch, remove the test cassette. Place the test cassette on a flat surface and unscrew the bottle of reagent and place it aside.

- 1 Take up the serum or plasma sample with the pipette.
- 2 Carefully put one (1) drop (20µl) of sample material into the sample well.
- 3 Add two (2) drops of the reagent from the bottle of reagent into the sample well. Ensure that no air bubbles are formed.

The liquid starts running up the test strip after a short time (< 90 seconds). If the fluid does not run up the test strips after 90 seconds, add an additional drop of the reagent into the sample well and/or **press with the tip of the pipette into the sample well to reactivate the run of the test.**

 After 15 minutes, the result can be read. **Do not read off results after 15 minutes! After this time, the test and its results are no longer considered valid.**

Whole blood

A blood sample should be used as soon as possible. Heparinised or EDTA blood may also be used.


The sample should be at room temperature and shaken well before use. Haemolysed or lipaemic samples are not recommended.

The use of whole blood samples may lead to a lower level of sensitivity. In case of a negative result with whole blood, if infection is suspected, the test should be repeated with a whole blood serum or plasma sample to obtain maximum sensitivity.

Test procedure

Open the pouch, remove the test cassette. Place the test cassette on a flat surface and unscrew the bottle of reagent and place it aside.

- 1 Take up the serum or plasma sample with the pipette.
- 2 Carefully put one (1) drop (20µl) of sample material into the sample well.
- 3 Add two (2) drops of the reagent from the bottle of reagent into the sample well. Ensure that no air bubbles are formed.


The liquid starts running up the test strip after a short time (< 90 seconds). If the fluid does not run up the test strips after 90 seconds, add an additional drop of the reagent into the sample well and/or **press with the tip of the pipette into the sample well to reactivate the run of the test.**

 After 15 minutes, the result can be read. **Do not read off results after 15 minutes! After this time, the test and its results are no longer considered valid.**

9. Test result

Positive result

The test results can be read after 15 minutes.


In case of a positive result, two red lines appear in the reaction field of the test cassette. A red line in the T-region (T) of the reaction field indicates a positive test result. Also a faint test line is considered as a positive test result. The second red line in the C-region (C) indicates the control line, which indicates the correct performance of the test.

The C-line is not a reference line and may have a different line intensity than the T-Line.

Note: a positive result is not conclusive but is an indicator in combination with a corresponding preliminary report. A

definitive diagnosis can only be established with the inclusion of additional clinical examinations.

Negative result

If only the Control Line C appears, it means a negative result. This indicates that no Leishmania, Ehrlichia or Anaplasma antibodies or dirofilaria antigens have been detected in the sample material, depending on the strip.

Note: Serological results must always be interpreted in combination with clinical results. Negative serum results do not rule out an infection because these results can be obtained at any stage of infection. Positive serum results in endemic areas may be due to a previous infection.

Invalid result

If **no** control line appears at C after the test has been performed, the test is invalid. In this case, it is likely that the test has not been performed correctly or that the expiry date has been exceeded. In this case, a new test must be performed.

If the problem persists, stop using the kit immediately and contact your dealer or our technical customer service department on **0034 958 41 27 76.**

10. Quality control

An internal control is included in the test cassette. The red line in the control region (C) is considered an internal control of the procedure. This line confirms that the sample volume was adequate, that the membrane was sufficiently soaked and that the procedure technique was correct.

11. Limitations

The test results should always be judged in connection with the anamnestic and clinical context. All test results obtained with whole blood samples, which do not correlate with the clinical picture should be repeated with a serum or plasma sample. All doubtful test results should be repeated after approximately 1 week or referred to the laboratory for certification. Positive results detected with the GaSS-Test FOuR VECt veterinary rapid tests must be confirmed with an additional analytical method.

There is a possibility that the test result may be distorted by technical errors, procedural faults or by substances or factors not mentioned here, but which influence the test. For questions, comments or technical queries, please contact our technical customer service department:

0034 958 41 27 76.

Grupo DAV SALUD
C/Tucumán nº6 Nave A-B,
18200 Maracena (Granada), Spain.

12. Performance characteristics

Test	Sensitivity	Specificity	Overall agreement	Comparative method
Leishmania Ac	97,78 %	100,00 %	98,36 %	IFI
Ehrlichia Ac	95,58 %	99,99 %	98,00 %	IFI
Anaplasma Ac	96,49 %	99,99 %	98,88 %	IFI
Dirofiladia Ag	96,72 %	99,99 %	98,62 %	ELISA

13. Literature

E Teske, F van Knapen, EGM Beijer, RJ Slappendel: " Risk of Infection with Leishmania spp. in the Canine Population in the Netherlands" in Acta Veterinaria Scandinavica 2002, 43:195-201.

Lindsay DS, Zajac AM, Barr SC: Leishmaniasis in American Foxhounds: An Emerging Zoonosis? Compend Cont Educ Pract Vet 24:304-312, 2002.

Torsten J. Naucke: „Leishmaniose, eine Tropenkrankheit und deren Vektoren (Diptera, Psychodidae, Phlebotominae) in Mitteleuropa“ in Denisia 0006, 2002 (163-178).

Al-Adhami B., Scandrett W.B., Lobanov V.A. & Gajadhar A.A. (2011): Serological cross-reactivity between Anaplasma marginale and an Ehrlichia species in naturally and experimentally infected cattle. Journal of Veterinary Diagnostic Investigation 23, 1181-1188.

Alleman, A. R. and Wamsley, H. L. : An update on anaplasmosis in dogs. *Vet Med* 2008, April: 212-220.

CDC: Ehrlichiosis and Anaplasmosis: 2008 Case Definition. US Department and Health and Human Services, CDC, National Notifiable Diseases Surveillance System (NNDSS): Atlanta, GA.

Dujardin J-C, Campino L, Carmen C, Dedet, J-P, Gradoni L, Soteriadou K, Mazeris A, Ozbil Y and Boela M. 2008. spread of vector-borne diseases and neglect of leishmaniasis, Europe. *Emerging Infectious Diseases*. *Emerging Infectious Diseases*. 14, no. 7.

Eckert, Johannes, Friedhoff, Karl Theodor, Zahner, Horst, Deplazes, Peter. Textbook of parasitology for veterinary medicine. 2nd, completely revised edition. Enke-Verlag, 2008
European Scientific Counsel Companion Animal Parasites: Control of vector-borne diseases in dogs and cats.

Maggi, R. G.; Birkenheuer, A. J.; Hegarty, B.C.; Bradley, J. M.; Levy, M. G.; Breitschwerdt, E. B. : Comparison of serological and molecular panels for diagnosis of vector-borne diseases in dogs. *Parasite Vectors*. 2014;7:127.

Pantchev, N.: Tick-borne travel infections in dogs: ehrlichiosis and babesiosis. In: *Veterinary Review*. Vol. 4, 2012, pp. 162-170.

Ranjbar-Bahadori, S. H.; Veshgini, A.; Shirani, D.; Eslami, A.; Mohieddin, H.; Shemshadi, B.; Masooleh, R. : Epidemiological aspects of canine dirofilariasis in the north of Iran. *Iran J Parasitol*. 2011 Mar;6(1):73-80.

Sainz, Á.; Roura, X.; Miró, G.; Estrada-Peña, A.; Kohn, B.; Harrus, S.; Solano-Gallego, L.: Guideline for veterinary practitioners on canine ehrlichiosis and anaplasmosis in Europe. *Parasite Vectors* 2015 Feb 4;8:75. doi: 10.1186/s13071-015-0649-0.

Self, S.; Liu, Y.; Nordone, S. K.; Yabsley, M. J.; Walden, H. S.; Lund, R. B.; Bowman, D. D.; Carpenter, C.; McMahan, C. S.; Gettings, J. R.: Canine vector-borne disease: mapping and the accuracy of forecasting using big data from the veterinary community. *Anim Health Res Rev*. 2019 Jun;20(1):47-60. doi: 10.1017/S1466252319000045.

Vieira, A. L.; João Vieira, M.; Oliveira, J. M.; Simões, A. R.; Diez-Baños, P.; Gestal, J. : Prevalence of canine heartworm (*Dirofilaria immitis*) disease in dogs of central Portugal. *Parasite*. 2014; 21:5. doi: 10.1051/parasite/2014003. epub 2014 Feb 19.

GASSET[®]
laboratorios

Grupo DAV SALUD

C/Tucumán nº6 Nave A-B, 18200 Maracena (Granada), España/Spain.

Tlf. Asistencia Técnica 958 41 27 76